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ELITE LC TRAY ROUND INSTRUCTIONS FOR USE - EN

Light-curing compound resin-based material for the preparation of individual trays (special trays), available in pre-formed plates. May be used in daylight or under artificial lighting for the amount of time required for modelling. Easy to use with high dimensional stability, rigidity and resistance to torsional stress. Plates are supplied in boxes to protect them against the light.

Leftover material may be stored in its box and used later.

Elite LC Tray ROUND may be polymerised with any kind of UV device (wavelength of 350-400 nm) or halogen device (wavelength of 420-480 nm).

Uses:

Individual trays (special trays)

Functional trays

Base plates for bite taking

Bases for assembly of teeth for dentures

Base plates for fixing registration devices

Fastening telescopic or secondary crowns in place

Instructions

Preparation

Trace the cutting proof for the individual spoon on the plaster model.

Discard any undercuts. Insert spacer blocks into the model if necessary, such as wax sheets. If no spacers are required, apply a thin layer of model insulation.

Application

Take the plate out of its box, closing the box again immediately.

Remove the protective film and put the resin plate in position by fitting it to the model. Cut it with a sharp instrument along the outline traced previously and model it. A handle for the tray may be made out of excess material.

Polymerisation

Put the model covered by the tray into a halogen or UV light polymerisation device. Hardening will take 3 to 5 minutes. Then remove the tray from the model, turn it over and expose it to the light for another 3 to 5 minutes. Repeated polymerisation will improve the material’s mechanical properties.

Important

A jelly-like sheen will form on the surface of the hardened resin after polymerisation due to inhibition caused by the presence of oxygen. It is important to remove this material with a tissue soaked in ethyl alcohol to permit effective adhesion.

Finishing

Remove any residues of wax with hot water.

Finish with resin cutters, smooth with abrasive paper and polish with a sander.

Caution: work with a dust suction system and be careful not to breathe in dust. Like all metacrylate resins, it may cause allergic reactions on sensitive skin; in the event of such a reaction, suspend use of the product.

Disinfection: compatible with disinfectants normally used on tray materials For optimum disinfection, use: Zeta 7 Solution or Zeta 7 spray (Zhermack).
Storage and duration
Elite LC Tray ROUND lasts 3 years if stored at room temperature and protected from light. Do not store at temperatures in excess of 25°C.

Disinfection: compatible with disinfectants normally used on tray materials For optimum disinfection, use: **Zeta 7 Solution** or **Zeta 7 spray** (Zhermack).

Storage and duration

Elite LC Tray ROUND lasts 3 years if stored at room temperature and protected from light. Do not store at temperatures in excess of 25°C.

IMPORTANT REMARKS: Consulting services on the use of our products, provided in any form, be that verbal, written or through demonstrations, are based upon the current state of dentistry and of our know-how. They are to be considered as non-binding information, also in relation to possible third party claims, and do not exonerate the user from personally checking whether the product is suitable for the planned use. Use and application by the user take place under no possible control by the Company, and are therefore under the user’s responsibility. Possible liability for damages is limited to the value of goods supplied by the Company and utilized by the user.

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<div></div>
IMPORTANT REMARKS: Consulting services on the use of our products, provided in any form, be that verbal, written or through demonstrations, are based upon the current state of dentistry and of our know-how. They are to be considered as non-binding information, also in relation to possible third party claims, and do not exonerate the user from personally checking whether the product is suitable for the planned use. Use and application by the user take place under no possible control by the Company, and are therefore under the user’s responsibility. Possible liability for damages is limited to the value of goods supplied by the Company and utilized by the user.
SOLO PER USO DENTALE
POUR L’USAGE DENTAIRE SEULEMENT

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