## ELITE VEST PLUS DE

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## ELITE VEST PLUS **EN**<sup>ⓘ</sup>

**Composition**

Phosphate bonded, carbon free precise universal investment for crown & bridge and pressable ceramic. Compatible with base- and precious metal alloys (except titanium) and pressable ceramic. Suitable for speed and overnight pre-heating technique.

Technical data	
- Powder/Liquid ratio (Elite Vest Liquid)	100 g / 24-26 ml
- Manual pre-mixing time	15 -30"
- Vacuum mixing time (320 rpm)	60"
- Vacuum maintaining time after mixing	15 -30"
- Processing time (20-23 °C)	6'
- Time for placing it in the furnace from mixing start	23 -25'
- Ideal usage temperature	20-23 °C
- Liquid storage temperature (Elite Vest Liquid)	> 5 °C
- Maximum pre-heating temperature	1200 °C

The indicated times may vary according to material and room temperature. High temperature reduces processing time and makes the material harden faster.

**Indications**

Suitable for ring and ring-less technique for casting crown and bridge made of base - and precious-metal alloys and pressable ceramics (except titanium).

Classification: UNI EN ISO 15912

Type 1, Class 1,2

**Mixing liquid**

Elite Vest Liquid 1000ml (universal for Elite Vest investment)

**Storage**

Store the product in a ventilated place away from light and at suggested temperature of:

- Elite Vest Plus 18-22°C (max 40°C)
- Elite Vest liquid 18-22°C (min 5°C)

Warning: Elite Vest liquid may freeze if exposed to temperatures below 2°C / 36°F. Once it freezes it can no longer be used.

<b>INSTRUCTIONS FOR USE</b>

**Using Elite Vest Plus with metal alloys**

**Dosage**

100 g Elite Vest Plus: 24-26ml Elite Vest Liquid

160 g Elite Vest Plus (bag) 39-41 ml Elite Vest Liquid

We recommend using a graduated measuring cup to dose Elite Vest liquid.

Note: for slow pre-heating, reduce the concentration of the liquid from 5 to 10% according to the alloy and application.

**Mixing ratio**

**Elite Vest liquid concentration recommended for a 160 g bag 40 ml liquid** (consider the data as a guideline).

		Wax crowns and bridges		Partial wax crowns, Inlays, Onlays		Telescopic resin crowns	
		Liquid Concentration		Liquid Concentration		Liquid Concentration	
High-gold alloys	Liquid	65%	26ml	50%	20ml	75%	30ml
	Distilled water		14ml		20ml		10ml
Low-gold alloys	Liquid	70%	28ml	60%	24ml	85%	34ml
	Distilled water		12ml		16ml		6ml
palladium alloys	Liquid	75%	30ml	65%	26ml	95%	38ml
	Distilled water		10ml		14ml		2ml
Base-metal alloys	Liquid	100%	40ml	90%	36ml	100%	40ml
	Distilled water		0ml		4ml		0ml

**General Elite Vest Liquid concentration indications per bag using metal alloys** (consider the data as a guideline).

% liquid concentration	160 g powder		320 g powder		480 g powder	
	Liquido ml	distilled water ml	Liquido ml	distilled water ml	Liquido ml	distilled water ml
100%	40	-	80	-	120	-
95%	38	38	38	2	76	4
90%	36	4	72	8	108	12
85%	34	6	68	12	102	18
80%	32	8	64	16	96	24
75%	30	10	60	20	90	30
70%	28	12	56	24	84	36
65%	26	14	52	28	78	42
60%	24	16	48	32	72	48
55%	22	18	44	36	66	54
50%	20	20	40	40	60	60

**Mixing**

Dose Elite Vest Plus and Elite Vest Liquid and pour the liquid and then the powder in a bowl; mix manually with a spatula for approximately 15-30". Then stir the mixture under vacuum for 1' using an automatic mixer (we recommend Zhermack Elite Mix). Once the mixture is complete, maintain the vacuum for other 15-30" before using the product. Note: ensure that the mixer creates the vacuum correctly.

**Processing**

Use metal or silicone cylinders. In the metal type we recommend inserting 1 expansion liner with vaseline (2 liners for X6 and X9 cylinders). Fill the cylinder pouring the Elite Vest Plus slowly to prevent bubbles. Use a low frequency vibrator to improve the smoothness of the investment.

Note: for X9 cylinders the speed process is not recommended

**Speed Pre-heating**

Once the mixture is set (23-25" after starting to mix) roughen the upper surface of the investment before pre-heating. Then place it in the pre-heated furnace at the final temperature (850 - 900 °C).

Note: gold alloys (Au) with pre-heating temperature ranging between 700 and 750 °C allow for better surface quality results by heating the cylinder investment at 850 - 900 °C first and then reducing the temperature of the furnace down to 700-750 °C.

Time at final temperature:

Cylinder	Position	The indicated times refer to one cylinder only. In the event that several cylinders are pre-heated at the same time, we recommend increasing the time by 10' for every additional cylinder.
X1	35'	
X3	45'	
X6	60'	

**Overnight Pre-heating**

For slow pre-heating, reduce the concentration of the liquid from 5 to 10% according to the alloy and application.

Once the mixture is set (25-30' after starting to mix) roughen the upper surface of the cylinder investment before pre-heating it. Then place the cylinder in the furnace at room temperature. Start the multi-stage program with the following

Position	Temperature	Holding time	Heat-rate
1	290°C	30-6	